



HAL
open science

Molecularly imprinted sol-gel polymers for the analysis of iprodione fungicide in wine: Synthesis in green solvent

Manal Bitar, Céline Lafarge, Nicolas Sok, Philippe Cayot, Elias Bou-Maroun

► To cite this version:

Manal Bitar, Céline Lafarge, Nicolas Sok, Philippe Cayot, Elias Bou-Maroun. Molecularly imprinted sol-gel polymers for the analysis of iprodione fungicide in wine: Synthesis in green solvent. *Food Chemistry*, 2019, 293, pp.226-232. 10.1016/j.foodchem.2019.04.108 . hal-02173340

HAL Id: hal-02173340

<https://u-bourgogne.hal.science/hal-02173340>

Submitted on 22 Oct 2021

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



Distributed under a Creative Commons Attribution - NonCommercial 4.0 International License

1 **Molecularly imprinted sol-gel polymers for the analysis of iprodione fungicide**
2 **in wine: synthesis in green solvent**

3 Manal BITAR^a, Céline LAFARGE^b, Nicolas SOK^{b,*}, Philippe CAYOT^b, Elias BOU-MAROUN^{b,*}

4
5 ^a *Université Libanaise, Faculté des Sciences IV, Laboratoire Energétique et Réactivité à l'Echelle*
6 *Nanométrique (EREN), Haouch El-Omara, Zahlé, Liban*

7 ^b *Univ. Bourgogne Franche-Comté, AgroSup Dijon, PAM UMR A 02.102, Procédés Alimentaires et*
8 *Microbiologiques, F-21000 Dijon, France*

9
10 **Abstract**

11 Iprodione is a fungicide widely used in viticulture in most agricultural countries. It was banned
12 recently in the European community because of its carcinogenic and endocrine disrupting characters.
13 In this work, a cheap analytical method able to monitor iprodione in a white wine was developed.
14 Molecularly imprinted sol-gel polymers (MIS) specific to iprodione and using green solvents were
15 synthesized. An experimental design having the following factors (solvent volume and crosslinker
16 quantity) was used to prepare an optimal MIS. In terms of selectivity, the optimal MIS showed the
17 best partition coefficient towards iprodione in a white wine containing four other competing
18 fungicides (procymidone, pyrimethanil, azoxystrobin and iprovalicarb). A solid phase extraction
19 method using the optimal MIS was optimized and applied to analyse iprodione in a white wine. Low
20 detection and quantification limits were reached 11.7 and 39.1 µg/L respectively.

21
22 **Keywords**

23 Molecularly imprinted silica, fungicide, wine, experimental design, solid phase extraction, cheap
24 method, green solvent

* Corresponding authors: elias.bou-maroun@agrosupdijon.fr, Phone: + 33 3 80 77 40 80 and
nicolas.sok@agrosupdijon.fr, Phone: + 33 3 80 77 40 48

25 **1. Introduction**

26 Iprodione is a dicarboximide fungicide that is widely used in agriculture around the world. It was
27 authorized to be used on vegetable crops, seed crops, fruit-tree crops, vines, ornamental crops,
28 grasses, fragrance, food, medicinal and culinary plants and field crops (ANSES, 2017).

29 According to a pesticide-monitoring programme performed by the French Ministry of Agriculture
30 (Cugier & Bruchet, 2005), iprodione was detected in 90% of the grapes analysed samples (47/52) and
31 in 100% of the wines analysed samples (43/43).

32 In 2008, iprodione was classified as carcinogen category 2 in accordance with (EC) regulation
33 N 1272/2008 of the European Parliament and of the council (Regulation EC 1272, 2008). Based on a
34 European Food Safety Authority study, the European Community concluded that there is a high risk
35 that the iprodione uses may expose groundwater to levels exceeding the drinking water limit of
36 0.1 µg/L (European Food Safety Authority (EFSA), 2016). An evaluation conducted by the French
37 Agency for Food, Environmental and Occupational Health and Safety (ANSES) in 2016 showed that
38 iprodione has endocrine disruptor properties. For all these reasons, the (EU) regulation N 2017/2091
39 decided not to renew the approval of iprodione. Any grace period accorded by the state members
40 expired on 5 June 2018 (Regulation EU 2091, 2017).

41 In this context, a reliable and cheap method of iprodione analysis should be used in routine in order
42 to check that agriculture products comply with this new regulation. Currently, chromatographic
43 methods coupled to mass spectrometry are used to analysis fungicides in food and environmental
44 samples (Castro et al., 2018; Elbashir & Aboul-Enein, 2018; Oliva et al., 2018). The big advantage of
45 these methods is their low detection limits. However, they are expensive and require high technical
46 experience for the lab staff. As an alternative, we proposed in a previous study a molecularly
47 imprinted solid phase extraction method coupled to high performance liquid chromatography with a
48 UV detection in order to analysis iprodione in wine (Bitar, Cayot, & Bou-Maroun, 2014). In this
49 previous work, an acrylate base molecularly imprinted polymer (MIP) was used.

50 Molecularly Imprinted Polymers are synthetic materials widely used in analytical chemistry as
51 preconcentration phase. They are shape memory materials. They act as a host for a guest target
52 thanks to a double chemical and steric complementarity. They are able to recognize the target in a
53 complex food (H. Li et al., 2018), **environmental matrices (Bakkour, Bolotin, Sellergren, & Hofstetter,**
54 **2018; Ruggieri et al., 2015) or biological samples (Gao et al., 2017; Wei, Mu, Huang, & Liu, 2017).**

55 Molecularly imprinted polymers are synthesized in the presence of a template. After solubilization
56 and interaction of a template molecule with a functional monomer in a solvent, a crosslinker is added
57 to form the three-dimensional structure of the polymer in the presence of an initiator. After washing
58 the polymer by an appropriate solvent, the template is removed leaving cavities highly specific to the
59 target.

60 Most of the imprinted polymers studied in the literature are acrylate based polymers thanks to the
61 high availability of commercial functional monomers and the ease of synthesis in organic solvents. In
62 this work, sol-gel molecularly imprinted polymers (MIS) are developed instead of acrylate based
63 molecularly imprinted polymer (MIP). As they are synthesized in a water/ethanol mixture, MIS
64 respect one of the principle of green chemistry by using safer solvents than those used in MIP
65 synthesis. Prat et al. studied solvents toxicity and elaborated a solvent classification from least toxic
66 to most toxic. Water and ethanol are among the 9 recommended solvents of the 51 studied (Prat,
67 Hayler, & Wells, 2014). Moreover, sol-gel MIS are more stable and specific towards the target species
68 than acrylate based MIP, also allowing faster diffusion of analytes (da Costa Silva & Augusto, 2006).

69 **Molecularly imprinted polymers were already used as solid phase extraction (SPE) material in order**
70 **to preconcentrate pesticides from real samples. Li et al. used MIS-SPE to extract carbendazim**
71 **fungicide from vegetables before its determination by UPLC-UV with a detection limit of 3 ng/L (S. Li,**
72 **Wu, Zhang, & Li, 2016). They used MIP-SPE to preconcentrate kelthane and pyridaben pesticides**
73 **from vegetables before they analysis by GC-MS with a detection limit of 1 and 3.2 pg/g respectively**
74 **(S. Li, Xu, Wu, & Luo, 2016).**

75 The first objective of this study was to obtain a sol-gel molecularly imprinted polymer (MIS) specific
76 to iprodione. In order to get an optimal MIS able to pre-concentrate iprodione, a 2-factors
77 experimental design was developed. In the MIS synthesis, the choice of the relative proportions of
78 the reagents is a crucial point. For this reason, the following 2-factors were studied: the solvent
79 volume and the crosslinker quantity. Selectivity was then checked using the optimal MIS in a white
80 wine containing iprodione and four other competing fungicides. A solid phase extraction (SPE)
81 method was developed using the optimal MIS to monitor iprodione in a white wine. The detection
82 and the quantification limits of this method were determined and compared to the previous method
83 based on MIP (Bitar et al., 2014) and to a SPE-LC-TOFMS method (Gilbert-López, García-Reyes,
84 Mezcuca, Molina-Díaz, & Fernández-Alba, 2007).

85

86 **2. Material and methods**

87 *2.1. Chemicals*

88 Iprodione (97%, CAS number 36734-19-7), Azoxystrobin (Analytical standard, CAS number 131860-
89 33-8), Iprovalicarb (Analytical standard, CAS number 140923-17-7), (3-Aminopropyl)trimethoxysilane
90 (APTMS 97%, CAS number 13822-56-5), tetraethoxysilane (TEOS \geq 99%, CAS number 78-10-4),
91 ammonium hydroxide (NH₄OH 28-30%, CAS number 1336-21-6), ethanol (\geq 99.8%, CAS number 64-
92 17-5), methanol (gradient grade, CAS number 67-56-1) and acetonitrile (\geq 99.9%, CAS number 75-05-
93 8) were purchased from Sigma Aldrich, France. Procymidone (\geq 98 %, CAS number 32809-16-8) and
94 pyrimethanil (\geq 98 %, CAS number 53112-28-0) were bought from Chemos GmbH, Regenstauf,
95 Germany. Water used in all experiments was deionised and obtained from an Elga Ionic system
96 PURELAB Option.

97

98 *2.2. Synthesis of the molecularly imprinted and non-imprinted sol-gel polymers*

99 The molecularly imprinted sol-gel polymers (MIS) were prepared at 40 °C in a thermostatic water
100 bath under magnetic stirring. The template molecule, iprodione was first solubilized in ethanol
101 (iprodione is sparingly soluble in water, 16 mg·L⁻¹ for pH≤8, and slightly soluble in ethanol, 25 g·L⁻¹;
102 logP = 3.1; pKa = 9.2, logD = 2.25 at pH=10). Then, water was added, followed by the monomer
103 APTMS (highly soluble in water, 1000 g·L⁻¹, and ethanol, 100 g·L⁻¹; pKa=10.6; from pH1 to 6, logD=-2;
104 logD=0.37 at pH10) and the crosslinker TEOS (TEOS is sparingly soluble in water, 0.4 g·L⁻¹, but reacts
105 with water and slowly hydrolyses; TEOS is miscible with ethanol; logP=3.81). Finally, NH₄OH was
106 introduced. The reaction mixture was left under stirring for 20 hours. The polymers were separated
107 from the liquid phase by centrifugation at 10 000 g for 10 minutes at room temperature. In order to
108 eliminate iprodione and free the specific cavities, the polymers were washed several times with
109 ethanol until iprodione was no longer detectable by reversed phase high performance liquid
110 chromatography (HPLC) in washing solvents. After washing, the polymers were dried for 6 h at 60 °C.
111 In parallel, NIS (Non Imprinted Silica) were synthesized under the same synthesis conditions as those
112 of MIS, but without using the template molecule. NIS served as control polymers.

113 In a typical synthesis, the following molar ratios were used. For 7.5 or 15 mL of solvent, 0.1 mmol of
114 iprodione, 0.4 mmol of APTMS, 2 or 6 mmol of TEOS and 10 mmol of NH₄OH. The proportion
115 ethanol/water was 90/10, v/v. In these conditions, the pH synthesis medium varies between 11.5
116 and 11.7. In the sol-gel process, the polymer formation pass through hydrolysis of alkoxy silane than
117 condensation. At low pH, the hydrolysis kinetic is fast and the condensation kinetic is slow. The
118 condensation starts when hydrolysis is completed. At high pH, condensation is faster than hydrolysis.
119 pH 11.5-11.7 was chosen to have a good compromise between the rate of hydrolysis and
120 condensation.

121 In **figure 1**, the different chemical reactions are presented at the different steps, at the MIS synthesis
122 step (A), washing step of the MIS (B) and and the mixing step with MIS and wine contaminated with
123 iprodione where the interactions between the MIS and iprodione are involved (C).

124

125 *2.3. Binding properties and imprinting factor of MIS/NIS towards iprodione*

126 10 mg of imprinted polymer (MIS) or non-imprinted polymer (NIS) were suspended in 20 mL (50/50,
127 v/v, ethanol/water) solutions containing 1 mmol·L⁻¹ of iprodione. The mixture was stirred for 1 h at
128 25 °C. After removal of the polymer particles by centrifugation at 10 000 g for 10 min at room
129 temperature, the supernatant was diluted 10 times, then filtered using a 2 µm filter and analysed by
130 HPLC for the determination of iprodione.

131 For HPLC iprodione analysis, a Shimadzu HPLC (LC-20AT) pump equipped with a UV-Visible detector
132 (SPD-20A), a Restek column (Pinacle II, C18 5 µm, 150 x 4.6 mm) and a pre-column (C18, 5 µm, 10 x
133 4 mm) were used. Separations were conducted at room temperature. Isocratic acetonitrile/water
134 (60/40, v/v) was used as mobile phase and the flow rate was 1mL/min. Iprodione detection was done
135 at 220 nm. A calibration curve was constructed over the range 0.1–500 µmol·L⁻¹ of iprodione. **Eight**
136 **standards were prepared in (50/50, v/v, ethanol/water) solutions, the calibration curve equation was**
137 **$y = 47820 x + 254129$ ($y = \text{peak area}$; $x = \text{concentration}$) and the correlation coefficient $R^2 = 0.9989$.**

138 The amount of iprodione retained by the polymer (B, mg of iprodione/g of polymer) was calculated
139 by difference between the initial iprodione concentration (C₀, mmol/L) and the free iprodione
140 concentration (F, mmol/L) in the supernatant.

141 The imprinting factor (IF) was calculated by the equation (Spivak, 2005):

142
$$IF = \frac{K(MIS)}{K(NIS)} = \frac{B(MIS) \times F(NIS)}{F(MIS) \times B(NIS)}$$

143

144 where K is the partition coefficient $K = B/F$. The imprinting factor evaluates the specific binding
145 property of each polymer towards iprodione. It takes into account the sorption of iprodione by a MIS

146 and its corresponding control polymer NIS. It allows to get rid of the non-specific interactions that
147 result from the iprodione adsorption on NIS.

148 Statistical analysis (ANOVA with Tukey's post hoc analysis) was done using (XLSTAT, Addinsoft)
149 software.

150

151 *2.4. Selectivity studies*

152 20 mg of MIS or NIS were suspended in 10 mL of a Burgundy white wine (Chardonnay, Charles
153 Renoir). The same wine was supplemented by the following fungicides: iprodione, procymidone,
154 pyrimethanil, azoxystrobin and iprovalicarb at 0.2 mmol·L⁻¹ of each fungicide (cf. figure S1 of the
155 supplementary material section for the fungicide chemical structures). 500 µL of the MIS or NIS
156 suspension was mixed with 500 µL of the supplemented wine in a 2 mL Eppendorf tube. The tubes
157 were then agitated for 3 h using a Stuart rotator SB2. The tubes were finally centrifuged for 10 min at
158 20 000 *g* at room temperature. The supernatant were analysed by HPLC using the same system
159 described in the paragraph 2.3 in order to determine the concentration of each of the studied
160 fungicides. Each experiment was repeated three times.

161 The following gradient was used to separate the 5 fungicides: (solvent A water, solvent B acetonitrile)
162 50% of B from 0 to 2 min, 70% of B from 2 to 8 min, 90% of B from 8 to 10 min, 50% of B from 10 to
163 12 min. Finally, a 2 min reconditioning step was applied at the initial conditions. A calibration curve
164 was constructed over the range 0.1–500 µmol·L⁻¹ of each fungicide. Eight standards were prepared in
165 (90/10, v/v, water/ethanol) solutions. The calibration curve equations and the correlation
166 coefficients were $y = 18648x + 18401$ ($R^2 = 0.9991$) for iprodione; $y = 10682x + 8412$ ($R^2 = 0.9993$) for
167 procymidone; $y = 7428x - 2736$ ($R^2 = 0.9996$) for pyrimethanil; $y = 30934x + 7228$ ($R^2 = 0.9994$) for
168 azoxystrobin and $y = 7217x - 3109$ ($R^2 = 0.9996$) for iprovalicarb.

169 Statistical analysis (ANOVA with Tukey's post hoc analysis) was done using (XLSTAT, Addinsoft)
170 software.

171

172 *2.5. Molecularly Imprinted Solid Phase Extraction (MISPE)*

173 Empty polypropylene SPE cartridge (1 mL, Supelco, Bellefonte, USA) was packed with 20 mg of MIS or
174 NIS previously suspended in acetonitrile. The polymers (particle size > 20 μm) were packed between
175 two polyethylene frits of 20 μm porosity at each end. The SPE cartridges were firstly preconditioned
176 with 5 mL of methanol or ethanol and then with 5 mL of (water/ethanol, v/v, 87/13). 6 mL white
177 wine sample (Chardonnay, Charles Renoir) spiked with ($5 \mu\text{mol}\cdot\text{L}^{-1} = 1.65 \mu\text{g}/\text{mL}$) of iprodione was
178 loaded. The cartridges were subsequently washed with 2 mL of washing solvents, mixture of
179 water/acetone at different proportions: 70/30, 50/50 or 25/75, v/v. The elution step was performed
180 using 1 mL of elution solvents: water/acetone, where the acetone proportion was varied between 75
181 and 100% (v/v). The effluents from the different SPE steps were analysed by HPLC in order to
182 quantify iprodione using the same conditions described in paragraph 2.3. A calibration curve was
183 constructed over the range $0.1\text{--}5 \mu\text{mol}\cdot\text{L}^{-1}$ of iprodione. All experiments were done in triplicate.

184

185 *2.6. Characterisation of the optimal sol-gel polymer*

186 Fourier transform infrared (FTIR) spectra were recorded on a Perkin Elmer spectrum 65 FT-IR
187 spectrometer in the range $4000\text{--}500 \text{ cm}^{-1}$ using attenuated total reflectance sampling. 64 scans with
188 a resolution of 4 cm^{-1} were applied.

189 For the surface morphological characterization and Energy-dispersive X-ray (EDX) spectrum, samples
190 were suspended in ethanol then a drop was placed on a silicon grid and examined in Hitachi Scanning
191 Electron Microscope (SEM) SU1510.

192 The size distribution of the optimal polymer was measured by laser diffraction (Malvern Mastersizer
193 3000) after suspending 20 mg of polymer powder in 10 mL of water and sonication for 2 min.

194

195 3. Results and discussion

196 3.1. Experimental design for the determination of the optimal polymer

197 The objective of the experimental design was to develop an optimal polymer for the
198 preconcentration of iprodione fungicide from wine samples. The optimal polymer must have the
199 highest adsorption capacity and the best imprinting factor toward iprodione. The following 2²
200 experimental design was applied (**table 1**).

201 **Table 1. Coded levels and corresponding factors of the 2² experimental design.**

Factors	X ₁	X ₂
	Solvent volume (V, mL)	Crosslinker quantity (Q, mmol)
Levels		
-1	7.5	2
1	15	6

202

203 The solvent volume (X₁ = V) and the crosslinker quantity (X₂ = Q) and were selected as the two factors
204 of the experimental design. Two levels coded (-1) and (+1) were considered for each factor. This
205 experimental design implied the synthesis of 4 MIS (MISa, MISb, MISc and MISd). In parallel to the
206 MIS preparation, 4 NIS were synthesized. They constituted control polymers for the corresponding
207 MIS.

208 The influence of X₁ and X₂ factors on the following two variables: the imprinting factor (Y₁ = IF) and
209 the iprodione adsorbed amount (Y₂ = B) was studied. Y₁ and Y₂ represent the responses of the
210 experimental design . The determination of IF and B was done three times for each polymer. The
211 results are compiled in **table 2**.

212

213 **Table 2. Experimental design factors (X) and responses (Y). Q (crosslinker quantity), S (solvent volume), IF**
 214 **(Imprinting Factor), B (adsorbed iprodione, mg of iprodione/g of polymer).**

Experiment N°	Polymer	X ₁ = V	X ₂ = Q	Y ₁ = IF	Y ₂ = B (mg/g)
1	MISa,1	-1	-1	3.07	15.03
2	MISa,2	-1	-1	3.37	13.74
3	MISa,3	-1	-1	3.22	14.38
4	MISb,1	-1	1	0.69	27.34
5	MISb,2	-1	1	0.84	26.85
6	MISb,3	-1	1	0.88	27.95
7	MISc,1	1	-1	1.58	12.49
8	MISc,2	1	-1	1.69	13.38
9	MIS,c3	1	-1	1.69	11.05
10	MISd,1	1	1	0.84	12.61
11	MISd,2	1	1	0.79	11.94
12	MISd,3	1	1	0.81	12.27

215

216 For each response, the calculated model from the experimental design is expressed as follows:

217 $Y (\text{response}) = b_0 + b_1.X_1 + b_2.X_2 + b_{12}.X_1X_2$

218 b_0 is the response average. b_1 and b_2 are coefficient of the respective factors X_1 and X_2 , they
 219 represent factors effect on Y. The equation resolution for the 2 responses Y_1 and Y_2 gave:

220 $Y_1 = 1.6253 - 0.3912 X_1 - 0.8154 X_2 + 1.5699 X_1X_2$

221 $Y_2 = 16.5901 - 4.2962 X_1 + 3.2405 X_2 - 13.0241 X_1X_2$

222 The significance limit for each response was calculated based on the standard deviation of all
 223 experiments "S", it was equal to 1.96S. The significance limit was respectively 0.188 and 1.620 for Y_1
 224 and Y_2 . For both responses, the interactions between the two factors (solvent volume and crosslinker
 225 quantity) were the most significant. **Results of ANOVA and regression applied to the response of the**
 226 **design of experiment are presented in tables S1 and S2 of the supplementary material section.**

227 **Figure 2** shows that as the solvent volume (V) increases during the synthesis step, the imprinting
 228 factor decreases or remain constant (fig.2A) and the quantity of adsorbed iprodione decreases
 229 (fig.2B). This observation could be explained by a dilution effect. The higher the solvent volume, the

230 lower the reactants meeting probability with monomers, the fewer the cavities specific to iprodione.
231 Therefore, there is a great interest in decreasing the solvent volume during the MIS synthesis.

232 Concerning the quantity of crosslinker (Q): When Q increases, the imprinting factor decreases (fig.2A)
233 but the iprodione binding property remain constant or increases (fig.2B). This observation could be
234 explained by the increasing of the non-specific interactions between the iprodione and the polymers.
235 In this case, the non-specific interactions are responsible for the amelioration of iprodione binding
236 properties. Therefore, there is no interest in increasing the crosslinker quantity during the MIP
237 synthesis.

238 **Table 2** shows that MISb (V = 7.5 mL and Q = 6 mmol) had the best iprodione binding property
239 (27.4 mg/g). However, it had the smallest imprinting factor (0.81). MISa (V = 7.5 mL and Q = 2 mmol)
240 had the best imprinting factor (3.2) and the second best iprodione binding property (14.4). Based on
241 the experimental design results, MISa was chosen as the optimal MIS.

242 The optimal MIS imprinting factor was 3.2. It is higher than the one we determined in a previous
243 study (IF = 2.4) (Bitar et al., 2015) that aimed to develop acrylate-based MIP specific to iprodione. On
244 the other hand, the optimal MIS had an adsorption capacity of 14.4 mg/g (mg of iprodione per g of
245 MIS). It is 33 times lower than that of the MIP, which was equal to 473 mg/g. This comparison of the
246 acrylate-based (MIP) and the silica-based (MIS) polymers shows that the MIS are more selective but
247 have a relatively low adsorption capacity which is in agreement with the following studies reports in
248 literature (Cummins, Duggan, & McLoughlin, 2005; Marx & Liron, 2001).

249

250 *3.2. Characterization of the optimal MISa*

251 The optimal MISa was characterized by Scanning Electron Microscopy, laser diffraction
252 particle size analyser and Fourier transform infrared (FTIR) spectroscopy. Results are
253 presented in supplementary materials (figures S2 to S5). SEM micrographs show that the

254 bigger particles size of MISa range from 18 to 61 μm . Particle size distribution of MISa show
255 that the average particle size is $13 \pm 3 \mu\text{m}$. The EDX spectrum show 3 major peaks
256 corresponding to oxygen (52.6%), silicon (29.9%) and carbon (16.6%). These percentages
257 show that (3-Aminopropyl)trimethoxysilane (APTMS) monomer reacts with
258 tetraethoxysilane (TEOS) to form the sol-gel polymer. FTIR spectrum of MISa show that all
259 the observed peaks are related to the absorption of the TEOS crosslinker and the APTMS
260 monomer after the sol-gel hydrolysis and condensation. The condensed APTMS have several
261 absorption bands: the ones at 2984 cm^{-1} and 790 cm^{-1} correspond to the C-H group and the
262 one at 960 cm^{-1} corresponds to the C-C group. The condensed TEOS has an absorption band
263 at 1060 cm^{-1} corresponding to the Si-O-Si group.

264

265 *3.3. Selectivity of the optimal polymer MISa*

266 The selectivity of the optimal polymer MISa was studied in a Burgundy white wine supplemented
267 with the five following competing fungicides: iprodione, procymidone, pyrimethanil, iprovalicarb and
268 azoxystrobin. The fungicides choice was based on a study done by the French Ministry of Agriculture
269 in 13 regions representing the main French vineyards. The chosen fungicides were detected in most
270 of the 1 222 wine analysed samples (Cugier & Bruchet, 2005).

271 The partition coefficient K (mL/g) was calculated for each fungicide. The fungicides partition was
272 determined between MISa or NISa and the white wine. K values were calculated using the equation
273 (Duan, Yi, Fang, Fan, & Wang, 2013):

$$274 \quad K = \frac{B}{F} \times \frac{V}{m}$$

275 where B is the fungicide amount retained by the polymer and F the free amount remaining in wine.
276 V (mL) is the wine volume and m (g) is the polymer mass. The obtained results are presented in
277 **figure 3.**

278 Iprodione had the highest K MIS value followed by procymidone. These values are coherent because
279 MISa was synthesized using iprodione as template. Procymidone has a close chemical structure to
280 iprodione compared to the other fungicides. The same order was maintained for K NIS values.
281 However, K NIS values were much lower than K MIS values. This shows that the molecular imprinting
282 process was successful.

283 The selectivity of MISa was assessed by the selectivity coefficient α defined as follows:

$$284 \quad \alpha = \frac{K_{MIS}(\text{fungicide})}{K_{MIS}(\text{template})}$$

285 Where the template is iprodione. The calculated α values are: 1.1, 1.2, 1.7 and 2.3. They correspond
286 respectively to procymidone, azoxystrobin, pyrimethanil and iprovalicarb. Procymidone, the
287 fungicide having the closest chemical structure to iprodione, had the lowest α value. Iprovalicarb,
288 with the most distant chemical structure from iprodione, had the highest α value.

289 The optimal MIS was able to interact with iprodione but to a lesser extent with other competing
290 fungicides in a white wine sample. In order to eliminate the non-specific interactions and to pre-
291 concentrate only iprodione, an optimisation of molecularly imprinted solid phase extraction was
292 performed.

293

294 **3.4. Determination of iprodione in white wine by Molecularly Imprinted Solid Phase**

295 *Extraction (MISPE) using the optimal MISa*

296 Before the application of MISa in the preconcentration of iprodione from a white wine using solid
297 phase extraction (SPE), the optimization of the solvents used during the different steps of the SPE

298 was done. MISa and NISa were packed into two different cartridges. The SPE cartridge filled with
 299 MISa was washed with 10 mL of acetonitrile in order to eliminate any residual iprodione. Methanol
 300 and ethanol was tested as pre-conditioning solvent. Methanol allowed better sorption of iprodione
 301 and gave a better sorption difference between MISa and NISa. The conditioning solvent was
 302 water/ethanol, 90/10, v/v. After loading the white wine supplemented with iprodione, three washing
 303 solvents were tested: water/acetone, v/v, 70/30, 50/50 and 25/75 respectively. The latter were
 304 chosen because it allowed the elution of the total amount of iprodione adsorbed on the NIS. Several
 305 elution solvents were used in order to elute all iprodione adsorbed on the MIS. For this reason, the
 306 acetone proportion was varied between 75% and 100%. Acetone 100% eliminated the total amount
 307 of iprodione adsorbed by MISa. The adsorbed and recovered amount of iprodione during the
 308 different SPE steps are compiled in **table 3**.

309 **Table 3. Adsorbed and recovered amount of iprodione from white wine using MISPE. Means \pm SD (n = 3).**

Polymer	Adsorbed (μg)	Washing recovery (μg)	Elution recovery (μg)	Total recovery (%)
MISa	8.55 \pm 0.85	2.95 \pm 0.21	4.43 \pm 0.82	86 \pm 10
NISa	3.71 \pm 0.73	3.11 \pm 0.41	0	84 \pm 9

310

311 MISa adsorbed twice as much iprodione as NISa. The washing step removed all of the non-specific
 312 interactions that correspond to the amount of iprodione adsorbed on NISa. The iprodione recovered
 313 by the elution step corresponds to the specific interactions. It represents 52% of the adsorbed
 314 iprodione. Acceptable total recovery percentages (higher than 80 %) were obtained for both MISa
 315 and NISa.

316 Blank white wine was analyzed by HPLC in order to check the specificity of the developed method.
 317 The blank chromatograms showed no interfering peaks at the retention time of iprodione. The limit
 318 of detection (LOD) and the limit of quantification (LOQ) of the MISPE method were calculated using

319 the following equations: $LOD = (3.3/6) (SD/S)$; $LOQ = (10/6) (SD/S)$ where SD is the standard
320 deviation of the HPLC response and S is the slope of the calibration curve. 6 corresponds to the
321 preconcentration factor of the MIS-SPE method. The calculated LOD and LOQ were respectively 11.7
322 and $39.1 \mu\text{g}\cdot\text{L}^{-1}$. The relative standard deviation of the method ($n = 3$) was 10 %. The quantification
323 limit of this method is lower than that determined by MIP-SPE-HPLC in our previous studies in white
324 wine ($LOQ = 422 \mu\text{g}\cdot\text{L}^{-1}$) (Bitar et al., 2014). It is comparable to that determined by SPE-LC-TOFMS in
325 fruit juices ($LOQ = 1.5 \mu\text{g}\cdot\text{L}^{-1}$) (Gilbert-López et al., 2007) but really less expensive.

326 Currently, the fungicides detection methods used are mainly based on mass spectrometric detection
327 and are relatively expensive (Gilbert-López et al., 2007; Hayward & Wong, 2009; Wong et al., 2010).
328 In parallel, ELISA methods using an immunoabsorbent are used (Koch, Stier, Senseman, Sobek, &
329 Kerns, 2013; Watanabe & Miyake, 2007). In addition to their high cost, these methods cannot be
330 used in extreme pH or temperature environments. Preconcentration methods using molecularly
331 imprinted polymers (MIPs) could replace the methods currently used. MIPs are inexpensive materials
332 and resist to extreme temperature and pH conditions (Haupt, 2001; Svenson & Nicholls, 2001).
333 The developed MIS-SPE-HPLC method has two major limitations: it is time consuming and cannot be
334 applied for fungicides that do not absorb in the UV-visible range.

335

336 4. Conclusion

337 Molecularly imprinted sol-gel silica (MIS) were developed in this work in order to conduct synthesis
338 in safe condition, and to preconcentrate iprodione fungicide from a Burgundy white wine using solid
339 phase extraction. The analysis was done by a cheap HPLC-UV method. To prepare the optimal MIS for
340 iprodione, an experimental design having the following 2 factors (solvent volume and crosslinker
341 quantity) was developed. The responses of the experimental design were the imprinting factor (IF)
342 and the binding capacity (BC). The optimal MIS had a higher IF but a lower BC than a previously

343 synthesized acrylate based polymer (MIP) specific to iprodione. The design of experiment showed
344 that the optimal MIS with the higher IF and BC, was obtained with low solvent volume and low
345 crosslinker quantity during its synthesis. The selectivity of the optimal MIS was checked in a white
346 wine containing iprodione and four other competing fungicides. The fungicide with the closest
347 chemical structure to iprodione gave the lowest selectivity factor and the fungicide with the most
348 distant chemical structure to iprodione gave the highest selectivity factor. The optimal MIS was used
349 as packing material in a solid phase extraction (SPE) to pre-concentrate iprodione from a white wine.
350 After the optimization of the different steps of the SPE, MISPE-HPLC-UV method was able to detect
351 iprodione in a white wine with a detection limit of $11.7 \mu\text{g}\cdot\text{L}^{-1}$. This method could be used to monitor
352 this forbidden fungicide in wine easily and with a cheap lab equipment such as an UV-VIS-HPLC.

353

354 **Acknowledgements**

355 This work was supported by the Regional Council of Burgundy and the European Union via FEDER
356 (Fonds Européen de Développement Régional) funding. We are grateful to AgroSup Dijon for co-
357 funding this research program. The authors thank Fleur-Anne RIO and Camille DETANG for their
358 contribution to the experimental work. **They thank Anne-Sophie VIAUX for the Laser diffraction
359 analysis and Laurence DUJOURDY for helpful discussion about the design of experiment.**

360

361 **References**

- 362 ANSES. (2017). Phytopharmacovigilance. Synthèse des données de surveillance. *French Agency for*
363 *Food, Environmental and Occupational Health & Safety*. Retrieved from
364 https://www.anses.fr/fr/system/files/Fiche_PPV_Iprodione.pdf
- 365 Bakkour, R., Bolotin, J., Sellergren, B., & Hofstetter, T. B. (2018). Molecularly Imprinted Polymers for
366 Compound-Specific Isotope Analysis of Polar Organic Micropollutants in Aquatic
367 Environments. *Analytical Chemistry*, *90*(12), 7292–7301.
368 <https://doi.org/10.1021/acs.analchem.8b00493>
- 369 Bitar, M., Cayot, P., & Bou-Maroun, E. (2014). Molecularly imprinted polymer solid phase extraction
370 of fungicides from wine samples. *Analytical Methods*, *6*(16), 6467–6472.
371 <https://doi.org/10.1039/C4AY00619D>
- 372 Bitar, M., Maalouly, J., Chebib, H., Lebrét, A., Cayot, P., & Bou-Maroun, E. (2015). Experimental
373 design approach in the synthesis of molecularly imprinted polymers specific for iprodione
374 fungicide. *Reactive and Functional Polymers*, *94*, 17–24.
375 <https://doi.org/10.1016/j.reactfunctpolym.2015.07.002>
- 376 Castro, G., Pérez-Mayán, L., Rodríguez-Cabo, T., Rodríguez, I., Ramil, M., & Cela, R. (2018).
377 Multianalyte, high-throughput liquid chromatography tandem mass spectrometry method
378 for the sensitive determination of fungicides and insecticides in wine. *Analytical and*
379 *Bioanalytical Chemistry*, *410*(3), 1139–1150. <https://doi.org/10.1007/s00216-017-0724-9>
- 380 Cugier, J. P., & Bruchet, S. (2005). Plan de surveillance résidus en viticulture. Campagnes viticoles
381 1990-2003. Direction générale de l'alimentation. Sous direction de la qualité et de la
382 protection des végétaux.
- 383 Cummins, W., Duggan, P., & McLoughlin, P. (2005). A comparative study of the potential of acrylic
384 and sol-gel polymers for molecular imprinting. *Analytica Chimica Acta*, *542*(1), 52–60.
385 <https://doi.org/10.1016/j.aca.2005.01.042>

386 da Costa Silva, R. G., & Augusto, F. (2006). Sol–gel molecular imprinted ormosil for solid-phase
387 extraction of methylxanthines. *Journal of Chromatography A*, *1114*(2), 216–223.
388 <https://doi.org/10.1016/j.chroma.2006.03.073>

389 Duan, Z., Yi, J., Fang, G., Fan, L., & Wang, S. (2013). A sensitive and selective imprinted solid phase
390 extraction coupled to HPLC for simultaneous detection of trace quinoxaline-2-carboxylic acid
391 and methyl-3-quinoxaline-2-carboxylic acid in animal muscles. *Food Chemistry*, *139*(1), 274–
392 280. <https://doi.org/10.1016/j.foodchem.2013.02.007>

393 Elbashir, A. A., & Aboul-Enein, H. Y. (2018). Application of gas and liquid chromatography coupled to
394 time-of-flight mass spectrometry in pesticides: Multiresidue analysis. *Biomedical*
395 *Chromatography*, *32*(2), e4038. <https://doi.org/10.1002/bmc.4038>

396 European Food Safety Authority (EFSA). (2016). Peer review of the pesticide risk assessment of the
397 active substance iprodione. *EFSA Journal*, *14*(11). <https://doi.org/10.2903/j.efsa.2016.4609>

398 Gao, D., Wang, D.-D., Zhang, Q., Yang, F.-Q., Xia, Z.-N., Zhang, Q.-H., & Yuan, C.-S. (2017). In Vivo
399 Selective Capture and Rapid Identification of Luteolin and Its Metabolites in Rat Livers by
400 Molecularly Imprinted Solid-Phase Microextraction. *Journal of Agricultural and Food*
401 *Chemistry*, *65*(6), 1158–1166. <https://doi.org/10.1021/acs.jafc.6b05269>

402 Gilbert-López, B., García-Reyes, J. F., Mezcua, M., Molina-Díaz, A., & Fernández-Alba, A. R. (2007).
403 Determination of Postharvest Fungicides in Fruit Juices by Solid-Phase Extraction Followed by
404 Liquid Chromatography Electrospray Time-of-Flight Mass Spectrometry. *Journal of*
405 *Agricultural and Food Chemistry*, *55*(26), 10548–10556. <https://doi.org/10.1021/jf072272s>

406 Haupt, K. (2001). Molecularly imprinted polymers in analytical chemistry. *Analyst*, *126*(6), 747–756.
407 <https://doi.org/10.1039/B102799A>

408 Hayward, D. G., & Wong, J. W. (2009). Organohalogen and Organophosphorous Pesticide Method for
409 Ginseng Root. A Comparison of Gas Chromatography-Single Quadrupole Mass Spectrometry
410 with High Resolution Time-of-Flight Mass Spectrometry. *Analytical Chemistry*, *81*(14), 5716–
411 5723.

412 Koch, P. L., Stier, J. C., Senseman, S. A., Sobek, S., & Kerns, J. P. (2013). Modification of a
413 commercially-available ELISA kit to determine chlorothalonil and iprodione concentration on
414 golf course turfgrass. *Crop Protection*, *54*, 35–42.
415 <https://doi.org/10.1016/j.cropro.2013.07.017>

416 Li, H., Zhao, L., Xu, Y., Zhou, T., Liu, H., Huang, N., ... Ding, L. (2018). Single-hole hollow molecularly
417 imprinted polymer embedded carbon dot for fast detection of tetracycline in honey. *Talanta*,
418 *185*, 542–549. <https://doi.org/10.1016/j.talanta.2018.04.024>

419 Li, S., Wu, X., Zhang, Q., & Li, P. (2016). Synergetic dual recognition and separation of the fungicide
420 carbendazim by using magnetic nanoparticles carrying a molecularly imprinted polymer and
421 immobilized β -cyclodextrin. *Microchimica Acta*, *183*(4), 1433–1439.
422 <https://doi.org/10.1007/s00604-016-1765-z>

423 Li, S., Xu, M., Wu, X., & Luo, J. (2016). Synergetic recognition and separation of kelthane and
424 pyridaben base on magnetic molecularly imprinted polymer nanospheres. *Journal of*
425 *Separation Science*, *39*(20), 4019–4026. <https://doi.org/10.1002/jssc.201600699>

426 Marx, S., & Liron, Z. (2001). Molecular Imprinting in Thin Films of Organic–Inorganic Hybrid Sol–Gel
427 and Acrylic Polymers. *Chemistry of Materials*, *13*(10), 3624–3630.
428 <https://doi.org/10.1021/cm000983z>

429 Oliva, J., Martínez, G., Cermeño, S., Motas, M., Barba, A., & Cámara, M. A. (2018). Influence of matrix
430 on the bioavailability of nine fungicides in wine grape and red wine. *European Food Research*
431 *and Technology*, *244*(6), 1083–1090. <https://doi.org/10.1007/s00217-017-3031-y>

432 Prat, D., Hayler, J., & Wells, A. (2014). A survey of solvent selection guides. *Green Chemistry*, *16*(10),
433 4546–4551. <https://doi.org/10.1039/C4GC01149J>

434 Regulation EC 1272. (2008). Regulation (EC) No 1272/2008 of the European Parliament and of the
435 Council of 16 December 2008 on classification, labelling and packaging of substances and
436 mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending

437 Regulation (EC) No 1907/2006 (Text with EEA relevance). *Official Journal of the European*
438 *Union, L 353/1, 31.12.2008*, pp. 1–1355.

439 Regulation EU 2091. (2017). Commission Implementing Regulation (EU) 2017/2091 of 14 November
440 2017 concerning the non-renewal of approval of the active substance iprodione, in
441 accordance with Regulation (EC) No 1107/2009 of the European Parliament and of the
442 Council concerning the placing of plant protection products on the market, and amending
443 Commission Implementing Regulation (EU) No 540/2011 (Text with EEA relevance.). *Official*
444 *Journal of the European Union, L 297/25, 15.11.2017*, pp. 25–27.

445 Ruggieri, F., D'Archivio, A. A., Camillo, D. D., Lozzi, L., Maggi, M. A., Mercurio, R., & Santucci, S.
446 (2015). Development of molecularly imprinted polymeric nanofibers by electrospinning and
447 applications to pesticide adsorption. *Journal of Separation Science, 38*(8), 1402–1410.
448 <https://doi.org/10.1002/jssc.201500033>

449 Spivak, D. A. (2005). Optimization, evaluation, and characterization of molecularly imprinted
450 polymers. *Advanced Drug Delivery Reviews, 57*(12), 1779–1794.
451 <https://doi.org/10.1016/j.addr.2005.07.012>

452 Svenson, J., & Nicholls, I. A. (2001). On the thermal and chemical stability of molecularly imprinted
453 polymers. *Analytica Chimica Acta, 435*(1), 19–24. [https://doi.org/10.1016/S0003-](https://doi.org/10.1016/S0003-2670(00)01396-9)
454 [2670\(00\)01396-9](https://doi.org/10.1016/S0003-2670(00)01396-9)

455 Watanabe, E., & Miyake, S. (2007). Immunoassay for iprodione: Key estimation for residue analysis
456 and method validation with chromatographic technique. *Analytica Chimica Acta, 583*(2),
457 370–376. <https://doi.org/10.1016/j.aca.2006.10.031>

458 Wei, Z.-H., Mu, L.-N., Huang, Y.-P., & Liu, Z.-S. (2017). Imprinted monoliths: Recent significant
459 progress in analysis field. *TrAC Trends in Analytical Chemistry, 86*, 84–92.
460 <https://doi.org/10.1016/j.trac.2016.10.009>

461 Wong, J. W., Zhang, K., Tech, K., Hayward, D. G., Krynitsky, A. J., Cassias, I., ... Brown, D. (2010).
462 Multiresidue Pesticide Analysis of Ginseng Powders Using Acetonitrile- or Acetone-Based

463 Extraction, Solid-Phase Extraction Cleanup, and Gas Chromatography–Mass
464 Spectrometry/Selective Ion Monitoring (GC-MS/SIM) or –Tandem Mass Spectrometry (GC-
465 MS/MS). *Journal of Agricultural and Food Chemistry*, 58(10), 5884–5896.
466 <https://doi.org/10.1021/jf903851h>

467
468

469 **Figure captions:**

470 **Figure 1.** Proposed chemical reactions involved in the MIS synthesis. A. Basic hydrolysis and poly-
471 condensation of tetraethoxysilane crosslinker and (3-Aminopropyl)trimethoxysilane functional
472 monomer in the presence of iprodione (in blue). B. Removal of iprodione from the polymer by
473 washing with ethanol. C. Rebinding of iprodione (in blue) in the specific cavities.

474

475 **Figure 2.** Imprinting factor \pm SD (n = 3) (A) and binding capacity \pm SD (n = 3) (B) of the MIPs towards
476 iprodione in accordance with of the experimental design. Q is the crosslinker quantity used during
477 the MIP synthesis. Statistical comparisons were done by ANOVA with Tukey's post hoc analysis.
478 Values significantly different (P<0.05) are indicated by different small letters.

479

480 **Figure 3.** Fungicides partition coefficients between MISa or NISa and a Burgundy white wine spiked
481 with iprodione, procymidone, azoxystrobin, pyrimethanil and iprovalicarb at 0.1 $\mu\text{mol/mL}$ of each
482 fungicide. Values are means \pm SD (n = 3/group) partition coefficient. Statistical comparisons were
483 done by ANOVA with Tukey's post hoc analysis. Values significantly different (P<0.05) between
484 fungicides are indicated by different letters. K MIS and K NIS values were processed separately.

485

486

Figure 1.




Figure 2.





Figure 3.

